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
Appendix 7: Determination of Cathodic Disbonding

A section of plastic tube with $d = 75 \text{ mm}$ and $h = 100 \text{ mm}$ is placed on the plate and impermeably bonded to the coated surface of each test specimen using an electrically non-conductive adhesive (silicone). The tube serves as a container for sodium chloride solution with a molar concentration of $c_{(\text{NaCl})} = 0.5 \text{ mol/l}$. Before bonding a fault with $d = 6 \text{ mm}$ is inserted in the area of the pipe into the coating and cutter.

arisation of the tube is achieved using a po- hereby the target potential is set at $U_{0, \text{NaCl}} =$ equivalent to $U_H = -1260 \text{ mV}$, $U_{\text{cathode total}} = -1500$

risation can also be achieved in a simplified a short circuit with a magnesium anode if the about $U_H = -1200 \text{ mV}$.

carried out for 30 days at a temperature of

In order to determine the depth of disbonding, six incisions are made radially starting at the fault in the coating, going down to the surface of the casting. One hour after completion of the test, the coating is lifted with a blade starting at the fault. The average depth of disbonding (= loss of adhesion of the coating) starting from the edge of the fault to the initial bond is then determined.

The following limiting values apply to the depth of disbonding:

Temperature	Test duration	Depth of disbonding
$23 \pm 2 \text{ °C}$	30 days	$\leq 10 \text{ mm}$

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Overview

Potentiostat

- Used for coating research experiments
- Manually operated
- Cathodic Disbondment Test
- Multichannel system
- Datalogger (optional)

Description

The **P-M 100** is a modular system. It is a pure potentiostat, galvanostat mode is not available.

The **P-M 100** is conscious a very simple device, which is used for series measurements. Here, the focus is on hold experiments (e.g. research of coatings). The settings have to be made manually. The measurement values can be read from the two digital displays for current and potential. The potentiostat has a maximum current of ± 100 mA with a compliance voltage of ± 11 V. It has two current ranges, namely 100 mA and 20 mA. A light-emitting diode is used as an indicator of the internal potential = OCP and as overload display. For the adjustment of potential a nominal voltage generator with ± 2 V is used.

P-M 100 is available as a single device in a case but also with up to 4 modules in one case. If needed, a data logger for data acquisition and storage is also available.

Accessories

Cathodic Disbondment Test Cell

Test cell for tests according to DIN Norm 30670, DIN EN ISO 15711, ISO 21809 Part 1 Annex H, ASTM G-8, ASTM G-42.